

WARRANTY REGISTRATION CARD

Dealer's portion
Name of product: _____
Name of owner: _____
Address: _____
Tel No. (Mobile) _____ (H) _____ (O) _____
Email: _____ Gender: M F
Dealer's Stamp: _____
Date of Purchase: / / year
Serial No.: _____

*IMPORTANT: To qualify for the warranty, please fill in this card and mail to us within 14 days from the date of purchase.


WARRANTY CERTIFICATE

2-YEAR WARRANTY

Owner's portion
Name of product: _____
Name of owner: _____
Address: _____
Tel No. (Mobile) _____ (H) _____ (O) _____
Email: _____ Gender: M F
Dealer's Stamp: _____
Date of Purchase: / / year
Serial No.: _____

*NOTE: Please produce this card together with the purchase receipt for warranty service.

■ APPEARANCE AND KEY FUNCTIONS OF THE THERMOMETER



■ INTRODUCTION

Thank you for choosing the FORA IR19 forehead_foot_skin_surface thermometer. Please read this instruction manual first so you can use this thermometer safely and correctly. Please keep this instruction manual for future reference. This innovative medical device uses advanced infrared (IR) technology to measure temperature instantly and accurately on the skin.

■ HOW DOES IT WORK

The FORA IR19 forehead_foot_skin_surface thermometer detects the heat generated by the skin over the temporal artery. FORA IR19 thermometer converts this heat into a temperature value in one second. The device is also suitable for surface temperature measurements such as feet, skin and other surfaces.

■ HIGHLIGHTED FEATURES

- Fast and accurate: meets the accuracy requirement specified in ASTM E1965-98.
- Non-contact measurement that is no probe cover is required, thus saving the cost of replacement.
- A probe that can be wiped using a clean cloth slightly dampened with alcohol or water.
- Automatically turn off if left unused over 30 seconds.
- Big button improves handling and activation.
- Memory function allows you to recall previous results up to 10 measurements.

■ IMPORTANT SAFETY INSTRUCTIONS

READ THIS BEFORE USING AND KEEP THESE INSTRUCTIONS

The following basic safety precautions should always be taken.

- Close supervision is necessary when the device is used by, on or near children, handicapped persons or invalids.
- Use the device only for the intended use described in this manual.
- Do not use the device if it is not working properly, or if it has suffered any damage.

■ CAUTIONS & WARNINGS

- As with any thermometer, proper technique is crucial to getting accurate temperature readings. Please read this manual thoroughly and carefully before using.
- Always store the thermometer in a cool and dry place: temperature between -20°C and 60°C (-4°F and 140°F); relative humidity less than 95%.
- Avoid direct sunlight.
- Avoid dropping the thermometer.
- Do not touch the lens.
- Do not disassemble the thermometer.
- Basic safety precautions should always be observed, especially when the thermometer is used on or near children or disabled persons.
- The use of this thermometer is not intended as a substitute for a consultation with your physician.

■ RESTRICTIONS OF USE

This thermometer is clinically proven to produce accurate temperature measurements. However, please be advised if you have the following situations: The accuracy cannot be ensured when the thermometer is not clean. Check that the probe is clean before taking a measurement.

■ USING THE DEVICE

Do not remove the front cover to measure.

■ Forehead measurement

Step 1
Press and release the ON/ MEMORY button. When ready, the thermometer beeps and displays the last measurement.




Step 2
Perform a forehead measurement with a non-contact distance within 3 cm. Press and hold the SCAN button. At the same time slowly slide the thermometer probe across the skin/surface (from one temple to another). As long as you are holding the scan button the thermometer will continue to measure. Therefore, once you release the button the thermometer will provide you with a result. To conduct a spot temperature measurement, simply press and immediately release the SCAN button. Spot measurements are possible for measurements on the feet, skin or measuring the temperature of other surfaces.

■ REPLACE THE BATTERY

The thermometer comes with one CR2032 lithium coin battery . Replace it when " " appears. Please follow the steps to replace a new battery.

Step 1
Push down the battery cover.(Figure 1)



Step 2
Place the new battery (with the "+" sign facing upwards) into the battery compartment and press it in until the battery is firmly secured. (Figure 2)

Step 3
Recrew the battery cover.

NOTE

- Although the thermometer works when " " appears, we still recommend that you change the battery to obtain an accurate result.
- Remove the battery if stored for a long period of time.
- The battery should be kept away out of child's reach. If they are swallowed, promptly see a doctor for help.

■ ABOUT NORMAL BODY TEMPERATURE & FEVER

The body temperature can vary with the individual. It also varies at different parts of the body and different time of the day. The following table shows the statistical normal ranges from different sites. Please keep in mind that temperatures measured from different sites, even at the same time, should not be directly compared.

Fever means that the body temperature is higher than normal. This symptom may be potentially caused by infection, overdressing or immunization. But some people may not undergo fever even when they are ill. These include, but are not limited to, infants not more than 3 months old, persons with compromised immune systems, persons taking antibiotics, steroids or antipyretics (aspirin, ibuprofen, acetaminophen), or persons with some chronic illnesses. Please consult your physician when you feel ill even if there is no fever.

Table 1-1 Body Site Normal Temperature Range

Site	Normal Temperature Range
Oral	0.6°C (1°F) or more above or below 37°C (98.6°F)
Rectal/ear	0.3°C to 0.6°C (0.5°F to 1°F) higher than oral temperature
Axillary (armpit)	0.3°C to 0.6°C (0.5°F to 1°F) lower than oral temperature

*1. Body Temperature at WebMD;
website: <http://firstaid.webmd.com/body-temperature> retrieved at 2010 Jan 7.

■ CARE & CLEANING

- The probe is waterproof. Please wipe it clean with a dry cloth or dampened with 70% alcohol.
- The body of the thermometer is not water-resistant. Please never put the thermometer under a running tap water or immerse it with any fluid. Use a soft and dry cloth to clean it. Do not use abrasive cleaners.
- Store the thermometer in a cool and dry location. Please keep it from dust and away from direct sunlight.

■ TROUBLESHOOTING

The table below shows problems you may encounter. Please follow "what to do" to resolve problems. If the problem still exists, please call your local dealer for help.

Message	What it means	What to do
	Room temperature is below 10°C.	Put the thermometer under operating temperature range of 10°C to 40°C (50°F to 104°F).
	Room temperature is above 40°C.	Put the thermometer under operating temperature range of 10°C to 40°C (50°F to 104°F).
	Error occurs in the meter.	Please put the thermometer under operating temperature range of 10°C to 40°C (50°F to 104°F). If the problem is unsolved, please contact the customer service line.
	Battery is low and " " is shown on LCD.	May continue to use the thermometer, but should change a new battery as soon as possible.
	Battery is too low to measure temperature correctly. No measurements will be taken.	Change a new battery immediately.

■ SYMBOL INFORMATION

Symbol	Referent	Symbol	Referent
	Consult instructions for use		CE mark
	Manufacturer		Type BF equipment
	Serial number		Do not use if package is damaged
	Caution, consult accompanying documents		

■ DISPLAY SCREEN

	Temperature reading
	Forehead temperature indicator
	Temperature falls within normal range
	Celsius scale / Fahrenheit Scale
	Memory mode
	Low battery power indicator
	High temperature indicator

■ SPECIFICATIONS

Model no.	FORA IR19
Dimensions	149 (L) x 44 (W) x 54.23 (H) mm
Weight	73.2g (with battery)
Displayed temperature range	22°C to 44°C (71.6°F to 111.2°F)
Display resolution	0.1°C / 0.1°F
Accuracy	±0.3°C (±0.5°F) for the range of 35.0°C (96.8°F) ~ 41.0°C (102.2°F)
Temperature unit	°C or °F
Operating temperature range	10°C to 40°C (50°F to 104°F)
Operating humidity	95% RH or less
Storage temperature range	-20°C to 60°C (-4°F to 140°F)
Storage humidity	95% RH or less
Memory capacity	10 measurements


ForaCare Suisse AG
CH-9000 St. Gallen
www.forecare.ch



311-1119100-015

FORA IR19

Forehead_Foot_Skin_Surface Thermometer



Owner's Booklet

Version 1.0 04,2013

■ RESULT INDICATORS

Message	What it means
	Temperature falls within measurement range: 22°C to 37.9°C (71.6°F to 100.2°F).
	Temperature falls within measurement range: 38°C to 44°C (100.4°F to 111.2°F).
	Temperature measurement falls outside the displayed temperature range: 22°C to 44°C (71.6°F to 111.2°F). Please follow this manual to take a reading again.

Tips for Measuring Temperature

- As with any type of thermometer, you may observe slight variations in consecutive measurements. It is recommended that you take 3 temperature readings and use the highest one for the following situations:
 - Infants not more than 3 months old.
 - Children not more than 3 years old who have a compromised immune system, where the presence / absence of fever is critical.
 - When you are learning to use the thermometer.
- Please wait at least 30 minutes after the following situations before taking a reading.
 - After exercising, swimming, or taking a bath. When exposed to extreme temperatures.
 - During eating and talking.

■ RECALLING THE MEMORY

Your thermometer stores the 10 most recent readings in the memory. Make sure the thermometer is in " OFF " mode before recall.



Step 2
Press the ON / MEMORY button to review the result. Each time you press the ON / MEMORY button, LCD will display stored temperature (latest result shown first), together with " " and number (from 1 to 10). When the memory is full, the oldest result is deleted as the new one added.

Step 3
The meter will automatically turn off if left without pressing any button within 30 seconds.

TERMS & CONDITIONS OF WARRANTY

- We warrant this product to be free of defects in workmanship and materials within the said warranty period on the warranty certificate.
- During the warranty period, if this product is found to be defective, you may bring it, together with the purchase receipt and Warranty Certificate, on a carry-in basis to our office during normal business hours for warranty service. We will then repair or replace defective parts or exchange the whole product as we may choose, at no charge to the original owner. After such repair, replacement or exchange, the product will be warranted for up to the remainder of the warranty period.
- This warranty is valid only if the Warranty Certificate and Warranty Registration Card are duly completed with date of purchase, serial number and dealer's stamp, and if the Warranty Registration Card is sent to our office not later than 14 days from the date of purchase.
- This warranty is void if this product has been repaired or serviced by unauthorized person. This warranty does not cover defects caused by misuse, abuse, accident, tampering, lack of reasonable care, fire or any other acts beyond human control.
- Except as stated in the above paragraphs, we disclaim all other warranties, implied or expressed, including the warranties of merchantability or fitness for a particular purpose with respect to the use of this product. We shall not be liable for any direct, consequential or incidental damages arising out of the use or inability to use this product.

PLEASE AFFIX STAMP

DEALER'S INFORMATION