

FORA[®]

Clinically
Validated

PO200

Fingertip Pulse Oximeter

 **Bluetooth[®]**




Comfortable-Clip
Design



One-button design,
easy to use



FEATURES



Model No.	FORA P0200
Dimension & Weight	H63 x W37x D32 mm, 40g without batteries
Connectivity	Bluetooth 4.0
Display	LCD
Battery Life	Batteries can be used continuously for 20 hours (for reference only, it depends on different brands of AAA alkaline batteries)
Power Source	Two 1.5V AAA alkaline batteries
Measurement Range	0% to 100%
Resolution	1%
Accuracy	100% ~ 80% ±2%; 79% ~ 70% ±3%; others are undefined.
Method	Dual wavelength LED

Pulse Rate

Measurement Range	30 to 250bpm
Resolution	1bpm
Accuracy	±1bpm or ±1%, whichever is greater
Operating Conditions	10° C to 40° C; below 95% R.H. (non-condensing)
Meter Storage / Transportation Conditions	-13° F to 158° F (-25° C to 70° C); Below 95% R.H. (non-condensing)

Classification

Type BF Applied part	Type BF Applied part
Safety	IEC60601-1
EMC	IEC60601-1-2
Harmonized Standard	ISO 80601-2-61:2011
Water-resistance	IP22
Mode of Operation	Spot Check / Monitoring

ForaCare Suisse AG

Neugasse 55, CH-9000 St. Gallen, Switzerland

TEL : +41-71-220-1001 FAX : +41-71-220-1075

URL : www.foracare.ch | foraconnect.ch

E-mail : chservice@foracare.com

FORA P0200_DM_Rev3.0_20201116



foracare.ch