

FORA[®] Safety Lancet

ACS037

 30G / 1.4mm

 28G / 1.4mm




- **Safe:** Single use, self-destructing, No possibility of reuse
- **Easy:** Intuitive 2-step operation. Blood sampling easier than ever before
- **Sterile:** Needles are sterilized with Gamma irradiation



FORA® Safety Lancet (ACS037) is designed to make blood sampling simple and easy. Touch activation means the lancet can be used in 3 easy steps.

TWIST



Twist off and discard lancet cap. DO NOT PULL.

PRESS







Press platform firmly against sample site. Lancet will automatically retract into body of device.

DISPOSE



Dispose of in a suitable sharps container.

Color	Gauge / Depth	Amount of bleeding
	30G / 1.4mm	
	28G / 1.4mm	



CE 0123 REF ACS037

FORA®
ForaCare Suisse AG

ForaCare Suisse AG
Neugasse 55, 9000 St. Gallen, Switzerland
TEL : +41-71-220-1001 FAX : +41-71-220-1075
URL : www.foracare.ch
E-mail : chservice@foracare.com

